

MVE TEC 3000

Quick Reference Guide

Multilingual

TEC 3000 Quick Reference Guide

Table of Contents

Product Identification	
Display / Control Panel	2
Bottom Panel / Electrical / Physical Connections	3
Dewar Plumbing Connections	4
Adjusting Temperature Alarm Settings	5
Adjusting Level and Level Alarm Settings	6-7
High Level Alarm Setting	
High Level Setting	
Low Level Setting	
Low Level Alarm Setting	
Adjusting Display and Output Settings	8
Calibration of Temperature Probes	8
Password and Security Setup	9-10
Alarms and Definitions	11
Contact Information	11
German table of contents	12
Italian table of contents	23
Spanish table of contents	34
French table of contents	45

NOTE: All MVE models are a Class 1, externally powered, continuous operation medical device. They are not suitable for use with flammable anesthetics. This equipment has been tested and found to comply with the limits for medical devices to IEC 601-1-2: [or EN 60601-1-1-2:2001 or Medical Device Directive 93/42/EEC].

NOTE: MVE liquid nitrogen freezers should be installed by an authorized MVE Distributor per the TEC 3000 Technical Manual, PN 13289499.

TEC 3000 Quick Reference Guide

Table 1: Front Panel Identification

1	Front Panel	The front panel is the user interface for the TEC 3000. All displays and controls are located on the front panel
2	Display	A 4 x 20 Liquid Crystal Display (LCD) shows the value of all current conditions including Temperature and Liquid Nitrogen Level. The display also shows any current alarm conditions that may exist.
3	Start Fill	Used to manually fill the freezer with liquid nitrogen. This button may also be pressed to clear “fog” from the storage area to increase visibility.
4	Stop Fill	Used to manually end the filling cycle. After pressing this button, the unit will not resume filling until the low level setting is reached, or the “START FILL” button is pressed.
5	Setup	Used to adjust the value of all user adjustable functions including Level, Level Alarms, Temperature Alarms, Passwords, Etc.
6	Alarm Mute	Used to silence the audible alarm. Also used to reset the latching alarm after the alarm condition is corrected.
7	Up Arrow	Used to increase number values during “Setup” modes. Press once to decrease incrementally. Hold the button to scroll through values quickly. May also be used to toggle “YES/NO” or “ON/OFF” values.
8	Down Arrow	Used to decrease number values during “Setup” modes. Press once to decrease incrementally. Hold the button to scroll through values quickly. May also be used to toggle “YES/NO” or “ON/OFF” values.
9	Escape	Used to exit any menu or setup function.
10	Enter	Used to Select any menu for editing or save any user adjustable setting.

TEC 3000 Quick Reference Guide



Table 2: Bottom Panel Identification

1	Temp A	Connection for Temperature A probe.
2	Temp B	Connection for Temperature B probe.
3	Power Connection	Connection for incoming power.
4	COM 1	Connection for RS-485 communications interface.
5	COM 2	Redundant Connection for RS-485 communications interface.
6	Manual Fill	Manual Fill button used to override the TEC 3000 software and manually fill the freezer.
7	Discreet Alarm Connections	These connections monitor the alarms for High Temp A, Low Level, Low Battery, and LN2 Supply. They can be wired to any alarm system that uses Normally Open (NO) contacts.
8	Global Alarm Connections	These connections monitor ALL alarms. They can be wired to any alarm system that uses Normally Open (NO) or Normally Closed (NC) contacts.
9	Main Harness Connection	Connection for the main wiring harness. All of the connections for the solenoid valves, lid switch, bypass sensor, etc. are located here.
10	Level Connection	Connection for level sensing system. Level sensing is achieved via differential pressure sensor, located inside the TEC 3000. The tube from the level sensing line on the freezer connects here. The connection is a 1/8" (3.2mm) OD hose barb.

Dewar Plumbing Connections

Connect a transfer line (included with freezer) from an LN₂ supply tank to the fill connection at the rear of the freezer. Optimum supply tank pressure is 22 to 35 psi (1.5 to 2.4 bar). Although the plumbing assembly has a 50 psi (3.45 bar) pressure relief device, it is recommended that the supply tank be pressurized below 35 psi (2.4 bar) to reduce the LN₂ “flash-off” rate during filling and to maximize the cryogenic valve life. The supply line can be insulated to minimize LN₂ transfer losses. After the transfer hose is securely coupled to the freezer and supply tank, ensure all connections are leak free by opening the valve of the LN₂ supply tank and apply a soap and water solution to each field joint. You should not see bubbles forming at any joint. Wipe away excess soap and water when finished. Before removing the transfer hose, ensure the LN₂ supply tank valve is closed. Slowly and carefully loosen the transfer hose connection to vent any remaining pressure in the line before disconnecting the hose.

TEC 3000 Quick Reference Guide

The following section describes how to adjust temperature alarm settings. At any time during the following procedure, the user may exit the level by pressing the “Escape” button repeatedly, until the display returns to the “monitor” display mode. After 30 seconds of inactivity, the controller will automatically return to the “monitor” display mode.

NOTE: Security Level 2 or higher is required to adjust temperature settings (see “Password and Security Setup” section for details).

To exit any menu screen and return to the previous menu, press “ESC” key.

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

High User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “ENTER”**

The display will read “Temperature Menus”.

Press Enter For
Temperature Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “Temp A Menu”.

NOTE: To adjust Temp B settings, press “Setup” instead of “Enter”.

Press Enter For
Temp A Menu
Or Press Setup For
Next Menu

4. **PRESS “ENTER”**

The display will read “Temp A Enabled”.

NOTE: Probe A may be disabled now by pressing the “▲ / ▼” keys. Press “Enter” to save settings.

Temp A
Enabled
Use ▲ / ▼ To Adjust
Press Enter To Save

5. **PRESS “SETUP”**

Until Display Reads “Temp A High Alarm”. Use the “▲ / ▼” keys to adjust the temperature. Hold either button to scroll quickly. Press “Enter” to save settings.

Temp A High Alarm
Use ▲ / ▼ To Adjust
Press Enter To Save
-110.0°C

6. **PRESS “SETUP”**

The display will read “Temp A Low Alarm”. Use the “▲ / ▼” keys to adjust the temperature. Hold either button to scroll quickly. Press “Enter” to save settings

Temp A Low Alarm
Use ▲ / ▼ To Adjust
Press Enter To Save
-200.0°C

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

TEC 3000 Quick Reference Guide

The following section describes how to adjust liquid nitrogen level settings and the high/low level alarms. NOTE: Security Level 2 or higher is required to adjust the Level and Level Alarm settings (See “Password and Security Setup” section for details).

NOTE: To exit any menu screen and return to the previous menu press “ESC” key.

HIGH LEVEL ALARM SETTING:

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move the cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Liquid Level Menus”.

Press Enter For
Liquid Level Menus
Or Press Setup For
Next Menu

2. **PRESS “ENTER”**

The display will read “High Level Alarm”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save the settings. NOTE: No level settings can be adjusted to less than .5 inches (14 mm) within each other.

High Level Alarm
8.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

HIGH LEVEL SETTING:

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Liquid Level Menus”.

Press Enter For
Liquid Level Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “High Level Alarm”.

High Level Alarm
8.0 In
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

The display will read “High Level Setpoint”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save settings.

High Level Setpoint
7.0 In
Use ▲ / ▼ To Adjust
Press Enter To Save

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

TEC 3000 Quick Reference Guide

Adjusting Level & Level Alarm Settings

NOTE: To exit any menu screen and return to the previous menu, press “ESC” key.

LOW LEVEL SETTING

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” arrow keys to scroll to the appropriate number. Press “Enter” to move the cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Liquid Level Menus”.

Press Enter For Liquid
Level Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “High Level Alarm”.

High Level Alarm
8.0 In
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

Until the display reads “Low Level Setpoint”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save settings.

Low Level Setpoint
5.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

LOW LEVEL ALARM SETTING:

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Liquid Level Menus”.

Press Enter For
Liquid Level Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “High Level Alarm”.

High Level Alarm
8.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

Until the display reads “Low Level Alarm”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save settings.

Low Level Alarm
4.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

TEC 3000 Quick Reference Guide

Adjusting Display And Output Settings

The units of measurement displayed by the TEC 3000 may be adjusted to accommodate the needs of the user. Temperature measurement may be displayed in degrees Kelvin (°K), degrees Celsius (°C), or degrees Fahrenheit (°F). The amount of liquid nitrogen in the freezer may be displayed in inches (in), millimeters (mm). In addition, the amount of liquid nitrogen consumed by the freezer (liquid usage) may be shown on the display.

NOTE: Security Level 1 is required to adjust the display and output settings (See “Password and Security Setup”). section for details).

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Display and Output”.

Press Enter For
Display And Output
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “Temperature Units”. Use the “▲ / ▼” keys to toggle between °C (Celsius), °F (Fahrenheit), or °K (Kelvin). Press “Enter” to save settings.

Temperature Units
°C
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

The display will read “Level Units”. Use the “▲ / ▼” keys to toggle between “in” (inches) or “mm” (millimeters). Press “Enter” to save settings.

Level Units
In
Use ▲ / ▼ To Adjust
Press Enter To Save

5. **PRESS “SETUP”**

The display will read “ Display Liquid Usage”. Use the “▲ / ▼” keys to select “ Yes” or “No”. Press “ Enter” to save settings.

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

Calibration of Temperature Probes

Each of the temperature sensor probes (A and B) used with the TEC 3000 have been calibrated at the factory using the “Low Temperature Range” method. This calibration method provides a level of accuracy of +/-1.8°F (+/-1°C) when operated in the altitude range of 1000ft. to 1500ft. (305m to 457m). Further calibration should not be required unless desired by the end user. Refer to the TEC 3000 Technical Manual for information on calibration methods and procedures.

TEC 3000 Quick Reference Guide

Password and Security Setup

The TEC 3000 can store up to 10 different passwords. Each password can be assigned its own security level ranging from Level 1 to Level 4. Table 3 below shows which settings can be changed with each security level. A security level of 4 is required to adjust any password. The default (or “Global”) password for the TEC 3000 is “3456”. All parameters may be adjusted by using this password. Record all passwords and security settings and store in a safe place. NOTE: MVE recommends changing the global password, as it is common to all units. If the global password has been forgotten, contact MVE Customer Service for details on how to reset passwords.

Table 3: Security Levels and Definitions

FEATURE	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4
Fill Start	X	X	X	X
Fill Stop	X	X	X	X
Alarm Mute	X	X	X	X
Change Display Units	X	X	X	X
Temp Settings		X	X	X
Level Settings		X	X	X
Time/Date		X	X	X
Calibration Probes		X	X	X
Change Languages		X	X	X
Hot Gas Bypass Settings		X	X	X
OFAF Setting			X	X
Communication Settings			X	X
Programming			X	X
Password Settings				X

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. If setting up passwords for the first time, use the default password of “3456”. Press “Enter” to move cursor to the next position in the password.

Higher Use Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “ENTER”**

The display will read “Temperature Menus”.

Press Enter For
Temperature Menus
Or Press Setup For
Next Menu

3. **PRESS “SETUP”**

Until the display reads “Password Menus”

Press Enter For
Password Menus
Or Press Setup For
Next Menu

4. **PRESS “ENTER”**

The display will read “Press Enter to change Global Password or press Setup for Next”.


Press Enter To
Change Global
Password Or Press
Setup For Next

NOTE: The default setting for the “Global” password is “3456”. The “Global” password may be used to change ALL settings for the TEC 3000.


TEC 3000 Quick Reference Guide

Password and Security Setup

5. **PRESS “SETUP”**
To scroll to the desired password number.
6. **PRESS “ENTER”**
To select the desired password number.
7. **PRESS “▲” or “▼”**
To scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.
8. **PRESS “ENTER”**
After the desired password is selected. The display will now read “Password X Level”. Use the “▲ / ▼” keys to adjust the security level.
9. **PRESS “ENTER”**
Display will read “Confirm New Password?”. Use the “▲ / ▼” keys to select “YES” or “NO”. Selecting “YES” will save the new password and security setting. Selecting “NO” will abort all changes.
10. **PRESS “ENTER”**
To save settings.

Password 1
Use ▲ / ▼ To Adjust
Press Enter For Next
XXXX

Password X Level
Use ▲ / ▼ To Adjust
Press Enter For Next
Level X

Confirm New
Password?
NO
Use ▲ / ▼ To Adjust

At this time, the selected password and security level has been set and saved. Follow steps 7 - 11 to set additional passwords and security levels, or press “Escape” repeatedly to return to the display screen. After 5 minutes of inactivity, the controller will automatically return to the “monitor” mode. Use Table 4 to record passwords for future reference.

TEC 3000 Quick Reference Guide

Alarms and Descriptions

Table 4: Alarms and Descriptions

Alarm Display	Description
High Temp A	The temperature of Probe A is above the user defined High Temperature setting.
High Temp B	The temperature of Probe B is above the user defined High Temperature setting.
Low Temp A	The temperature of Probe A is below the user defined Low Temperature setting.
Low Temp B	The temperature of Probe B is below the user defined Low Temperature setting.
High Level	The depth of LN2 inside the freezer is above the user defined High level setting.
Low Level	The depth of LN2 inside the freezer is below the user defined Low level setting.
Usage Warning	The consumption of LN2 has doubled (See TEC 3000 Technical Reference Manual for details).
Usage Alarm	The consumption of LN2 has increased by a factor of 5 (See TEC 3000 Technical Reference Manual for details).
Fill Time	The amount of time required to complete a fill cycle exceeds the user defined Fill Time setting.
Bypass Time	The amount of time required to complete a bypass cycle exceeds the user defined Bypass Time setting.
Temp A Calibration	The temperature of Probe A is lower than absolute zero.
Temp B Calibration	The temperature of Probe B is lower than absolute zero.
Bypass Calibration	The temperature of the Bypass Probe is lower than absolute zero.
Low Battery	The voltage of the back up batteries has dropped below 21 volts.
Power Failure	The primary power has been disconnected for at least 60 minutes.
Lid Open	The lid on the freezer has been open longer than the user specified time.
Communication Loss	The controller has lost communications with the display.

If any alarms occur, contact your authorized MVE Distributor or customer / technical service.

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Ball Ground, GA 30107, U.S.A.