



MVE 1879P-150  
MVE 1879P-190  
MVE 1892P-150  
SQUARE RACK LAYOUT




STANDARD



OPTIONAL


54 LARGE RACKS  
30 MINI RACKS  
13-2 RACKS EXCEPT 1892P-150  
WHICH HOLDS 15-2 RACKS

MVE 1879P-190  
BAXTER 4R9951 LAYOUT




492 FRAMES  
2952 CASSETTES

MVE 1879P-190  
BAXTER 4R9953 LAYOUT




396 FRAMES  
1584 CASSETTES

MVE 1879P-190  
BAXTER 4R9955 LAYOUT



276 FRAMES  
1104 CASSETTES


MVE 1879P-190  
GAMBRO DF200 LAYOUT



240 FRAMES  
960 CASSETTES




MVE 1879P-190  
GAMBRO DF700 LAYOUT



126 FRAMES  
504 CASSETTES

MVE 1879P-190  
MEDSEP LAYOUT



838 FRAMES  
5866 CASSETTES