



MVE 815P-150
MVE 815P-190
MVE 818P-150
MVE 819P-150
SQUARE RACK LAYOUT




12 LARGE RACKS
4 MINI RACKS
12-2 RACKS IN 815
14-2 RACKS IN 818
15-2 RACKS IN 819
15,600 VIALS IN 815
18,200 VIALS IN 818
19,500 VIALS IN 819

MVE 819P-150
BAXTER 4R9951 LAYOUT




CAPACITY
128 FRAMES
1,024 CANISTERS

MVE 819P-150
BAXTER 4R9953 LAYOUT




CAPACITY
104 FRAMES
520 CANISTERS

MVE 819P-150
BAXTER 4R9955 LAYOUT




CAPACITY
76 FRAMES
380 CANISTERS

MVE 819P-150
GAMBRO DF200 LAYOUT



CAPACITY
59 FRAMES
295 CANISTERS

MVE 818P-150
MEDSEP LAYOUT



CAPACITY
204 FRAMES
1,836 CANISTERS