



MVE 1536P-150
MVE 1536P-190
MVE 1539P-150
SQUARE RACK LAYOUT




24 LARGE RACKS
16 MINI RACKS
13-2 IN 1536
14-2 IN 1539
36,400 VIALS IN 1536
39,200 VIALS IN 1539

MVE 1539P-150
BAXTER 4R9951 LAYOUT




248 FRAMES
1736 CASSETTES

MVE 1539P-150
BAXTER 4R9953 LAYOUT




203 FRAMES
812 CASSETTES

MVE 1539P-150
BAXTER 4R9955 LAYOUT




152 FRAMES
608 CASSETTES

MVE 1539P-150
GAMBRO DF200 LAYOUT




124 FRAMES
496 CASSETTES

MVE 1539P-150
GAMBRO DF700 LAYOUT



64 FRAMES
256 CASSETTES

MVE 1539P-150
MEDSEP LAYOUT



383 FRAMES
3064 CASSETTES