



MVE 1894R-190  
BAXTER 4R9951 LAYOUT




490, 8 PLACE FRAMES  
3920 CANISTERS

MVE 1894R-190  
BAXTER 4R9953 LAYOUT




402, 5 PLACE FRAMES  
2010 CANISTERS

MVE 1894R-190  
BAXTER 4R9955 LAYOUT




310, 5 PLACE FRAMES  
1550 CANISTERS

MVE 1894R-190  
GAMBRO DF200 LAYOUT




246 FRAMES  
1230 CANISTERS

MVE 1894R-190  
GAMBRO DF700 LAYOUT



136 FRAMES  
680 CANISTERS


MVE 1894R-190  
MEDSEP LAYOUT



804, 8 PLACE FRAMES  
6432 CANISTERS

MVE 1881R-150  
MVE 1894R-190  
SQUARE RACK LAYOUT

Ø25.000"  
NECK ID



60 LARGE RACKS TOTAL  
12 MINI RACKS TOTAL  
13-2 RACKS IN 1881  
15-2 RACKS IN 1894  
81,900 VIALS IN 1881  
94,500 VIALS IN 1894