



MVE 616  
BAXTER 4R9951 LAYOUT




1 DIVIDER  
102 FRAMES  
612 CASSETTES

MVE 616  
BAXTER 4R9953 LAYOUT




1 DIVIDER  
84 FRAMES  
336 CASSETTES

MVE 616  
BAXTER 4R9955 LAYOUT




1 DIVIDER  
62 FRAMES  
248 CASSETTES

MVE 616  
GAMBRO DF200 LAYOUT




1 DIVIDER  
50 FRAMES  
200 CASSETTES

MVE 616  
GAMBRO DF700 LAYOUT





1 DIVIDER  
29 FRAMES  
116 CASSETTES

# MVE 616 MEDSEP LAYOUT




1 DIVIDER P/N 11742704  
196 FRAMES  
1372 CASSETTES

# MVE 616 SQUARE RACK LAYOUT



STANDARD



OPTIONAL

CAPACITY  
12, 13-2 LARGE RACKS 156, 100 CELL BOXES  
4, 13-2 MINI RACKS 52, 25 CELL BOXES  
16,900 VIALS