



MVE 816P-2T-190
SUC-1 LAYOUT



TOP TRAY
51 BOXES

1,173 C-1 CANES, 7,038 1.0ml VIALS
1,224 C-2 CANES, 7,344 1.2ml VIALS
816 C-7 CANES, 4,080 2.0ml VIALS
15,810 1/2cc STRAWS



BOTTOM TRAY
64 BOXES

1,472 C-1 CANES, 8,832 1.0ml VIALS
1,536 C-2 CANES, 9,216 1.2ml VIALS
1,024 C-7 CANES, 5,120 2.0ml VIALS
19,840 1/2cc STRAWS

TOTAL

115 BOXES

2,645 C-1 CANES, 15,870 1.0ml VIALS
2,760 C-2 CANES, 16,560 1.2ml VIALS
1,840 C-7 CANES, 9,200 2.0ml VIALS
35,650 1/2cc STRAWS