



MVE 205
SQUARE RACK LAYOUT




13-2 RACKS
4 LARGE (52 BOXES)
5200 VIALS

MVE 205
BAXTER 4R9951 LAYOUT




CAPACITY
1 DIVIDER
38 FRAMES
228 CASSETTES

MVE 205
BAXTER 4R9953 LAYOUT




32 FRAMES
128 CANISTERS

MVE 205
BAXTER 4R9955 LAYOUT




24 FRAMES
96 CANISTERS

MVE 2205
GAMBRO DF200 LAYOUT




20 FRAMES
80 CASSETTES

MVE 205
GAMBRO DF700 LAYOUT



10 FRAMES
40 CASSETTES

MVE 205 MEDSEP CANISTER LAYOUT



P/N	DESCRIPTION	QTY
11794176	DIVIDER	1
11044623	FRAMES	77
11856322	CANISTERS	539